



iWHITE
instant teeth whitening

Charcoal-infused Sonic Toothbrush for superior whitening



❖ Sonic Toothbrush

Product

The iWhite Sonic Toothbrush with charcoal bristles sweeps away plaque and gives you an amazing, beautiful, white smile.

The 40 000 high-speed sonic vibrations whip up your toothpaste into plaque fighting bubbles. These bubbles will drive deep between your teeth and along the gum line for exceptional cleaning. When used in combination with the **iWhite Dark Stains toothpaste**, it removes even better plaque and superficial stains from coffee, tea, wine and tobacco.

This electric toothbrush allows you to customize your brushing with a choice of 5 different modes and contains a rechargeable battery with 3 months of battery life.

Features



Charcoal

The unique charcoal bristles help to remove stains on the teeth to give you a beautiful smile.



Superior vibrations

The iWhite Sonic charcoal-infused toothbrush vibrates up to 40,000 times per minute to remove more dental plaque.



Timer

It only takes 2 minutes to clean your teeth thoroughly. The timer gives a signal each 30 seconds to change the quadrant. After 2 minutes, the toothbrush stops automatically.



5 brushing functions

1. Clean: the ultimate plaque removal
2. Sensitive: extra gentle 2-minute mode for sensitive teeth and gums
3. Gum Care: long amplitude vibration for gum massage
4. White: sonic vibration with alternate amplitudes to remove plaque and stains and whiten teeth
5. Polish: strong 2-minute vibration for teeth whitening and polish




3-month battery life

USB charging toothbrush, convenient to use at home or while travelling. Can be charged by PC, power bank, car or phone charger. If fully charged the battery will last up to 3 months (if used 4 minutes a day).

Instructions

1. Choose your favorite brushing mode.
2. Wet the bristles and apply a small amount of **iWhite Dark Stains toothpaste** on your toothbrush.
3. Place the toothbrush bristles against the teeth (do not press hard). To get a better cleaning, brush your teeth at a slight angle (45 degrees).
4. Gently keep the bristles placed on the teeth and in the gumline. Brush your teeth with small back and forth motion.
5. After you have completed brushing, the brush will turn off automatically after 2 minutes. If you would like to stop earlier, press the power button for 2 seconds.



Warning

Children aged below 12 years old are not recommended to use this product. Replace iWhite charcoal-infused brush heads every 3 months to achieve optimal results. Only iWhite replacement brush heads can be used on this device. Read carefully the internal leaflet for complete instructions for use, precautions and further information.

Content

1 Electrical toothbrush | 1 USB charger | 1 Replaceable head

Use the iWhite Sonic Toothbrush in combination with our other iWhite products for long-lasting whitening results.

Sylphar nv | Xavier De Cocklaan 42 | B-9831 Deurle, Belgium
Tel +32 9 321 70 80 | Fax +32 9 321 70 90 | www.iwhiteinstant.com

Sylphar